


Anti-Calponin-1 Rabbit Polyclonal Antibody For Wb, Elisa - P51911

Item Number: CM0018808



Introduction

Rabbit-derived polyclonal antibody targeting Calponin-1, raised against recombinant fusion protein. Ideal for Western blot and ELISA applications in smooth muscle research. Reacts with human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	Calponin Rabbit pAb
Remarks/Alias	SMCC; Sm-Calp; HEL-S-14; Calponin
Species	Human
Gene ID (Human)	1264
Gene ID	1264
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 200-297 of human Calponin (NP_001290.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51911
Protein Weight	33kDa
Shipping	Ice bag