


Anti-Calreticulin Monoclonal Antibody [Ko Validated] For Wb, Ihc-P, If/Icc, Elisa - P27797

Item Number: CM0008871



Introduction

Anti-Calreticulin mouse monoclonal antibody, KO validated, for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets CALR, a 47kDa ER chaperone involved in protein folding and calcium homeostasis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Calreticulin Mouse mAb
Remarks/Alias	RO; CRT; SSA; cC1qR; HEL-S-99n; Calreticulin
Species	Human
GeneID (Human)	811
GeneID	811
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-203 of human Calreticulin (NP_004334.1)
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27797
Protein Weight	47kDa
Shipping	Ice bag