


Anti-Calreticulin Rabbit Monoclonal Antibody (Ko Validated) For Wb, Ihc-P, Elisa - P27797

Item Number: CM0002119



Introduction

KO-validated rabbit monoclonal antibody targeting human Calreticulin (CRT/Erp60). Reacts with human, mouse, and rat. Confirmed for WB, IHC-P, and ELISA. ER chaperone and calcium-binding protein involved in glycoprotein folding.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Calreticulin Rabbit mAb
Remarks/Alias	RO; CRT; SSA; cC1qR; HEL-S-99n; [KD Validated] Calreticulin
Species	Human
Gene ID (Human)	811
Gene ID	811
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-203 of human Calreticulin (NP_004334.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27797
Protein Weight	47kDa
Shipping	Ice bag