


Anti-Calreticulin Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P27797

Item Number: CM0017778



Introduction

Calreticulin rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat calreticulin. Ideal for studying ER chaperone function, calcium homeostasis, and protein folding.

[Learn More](#)

Parameter	Value
Product Name	Calreticulin Rabbit pAb
Remarks/Alias	RO; CRT; SSA; cC1qR; HEL-S-99n; Calreticulin
Species	Human
GeneID (human)	811
GeneID	811
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-203 of human Calreticulin (NP_004334.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27797
Protein Weight	47kDa
Shipping	Ice bag