

Anti-Calretinin Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P22676

Item Number: CM0004135

Introduction

Rabbit monoclonal antibody targeting calretinin (CALB2), a 32 kDa calcium-binding protein. Validated in WB, IF-P, IHC-P, ELISA; cross-reacts with human, mouse, rat. Ideal for neuroscience and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Calretinin Rabbit mAb
Remarks/Alias	CR; CAL2; CAB29; Calretinin
Species	Human
GenelD (human)	794
GenelD	794
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-271 of human Calretinin (NP_001731.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22676
Protein Weight	32kDa
Shipping	Ice bag