

Anti-Calsequestrin 1 Monoclonal Antibody For Wb, If-P, Elisa - P31415

Item Number: CM0002789

Introduction

Rabbit monoclonal antibody against human Calsequestrin 1 (CASQ1), suitable for WB, IF-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying calcium storage and muscle contraction.

[Learn More](#)

Parameter	Value
Product Name	Calsequestrin 1 Rabbit mAb
Remarks/Alias	CASQ; CSQ1; PDIB1; VMCQA; Calsequestrin 1
Species	Human
Gene ID (Human)	844
Gene ID	844
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human Calsequestrin 1 (NP_001222.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P31415
Protein Weight	45kDa
Shipping	Ice bag