

Anti-Calu Rabbit Polyclonal Antibody For If/Icc, Elisa - O43852

Item Number: CM0016544

Introduction

Calumenin (CALU) rabbit polyclonal antibody for IF/ICC and ELISA, cross-reacts with human and rat. Ideal for studying vitamin K-dependent carboxylation and calcium binding.

[Learn More](#)

Parameter	Value
Product Name	CALU Rabbit pAb
Remarks/Alias	CALU; calumenin
Species	Human
GeneID (Human)	813
GeneID	813
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-315 of human CALU (NP_001210.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	O43852
Protein Weight	37kDa
Shipping	Ice bag