

Anti-Calumenin Monoclonal Antibody For Wb, Ihc-P, Elisa - O43852

Item Number: CM0009027

Introduction

High-quality rabbit monoclonal antibody against human Calumenin (CALU), validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for vitamin K-dependent carboxylation and calcium-binding studies.

[Learn More](#)

Parameter	Value
Product Name	Calumenin Rabbit mAb
Remarks/Alias	CALU; calumenin; Calumenin
Species	Human
Gene ID (Human)	813
Gene ID	813
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Calumenin (O43852).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43852
Protein Weight	37kDa
Shipping	Ice bag