


Anti-Camk1D Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Iu85

Item Number: CM0030608



Introduction

High-quality CAMK1D Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with Human, Mouse, Rat. Expected molecular weight 43kDa. Ideal for calcium signaling, kinase function, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	CAMK1D Rabbit pAb
Remarks/Alias	CKLIK; CaM-K1; CaMKID; CAMK1D
Species	Human
GeneID (Human)	57118
GeneID	57118
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human CAMK1D (NP_705718.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IU85
Protein Weight	43kDa
Shipping	Ice bag