


Anti-Camk1D Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Iu85

Item Number: CM0013977



Introduction

Rabbit polyclonal antibody against human CAMK1D, validated in WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for Ca²⁺/CaM kinase signaling and neutrophil function research.

[Learn More](#)

Parameter	Value
Product Name	CAMK1D Rabbit pAb
Remarks/Alias	CKLiK; CaM-K1; CaMKiD; CAMK1D
Species	Human
Gene ID (Human)	57118
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human CAMK1D (NP_705718.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IU85
Protein Weight	43kDa
Shipping	Ice bag