

Anti-Camk1G Polyclonal Antibody For Wb, Elisa - Q96Nx5

Item Number: CM0025173



Introduction

Rabbit polyclonal antibody targeting human CAMK1G (CaM kinase I gamma), validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studying calcium/calmodulin-dependent signaling and CREB1 phosphorylation in neuronal and peripheral tissues.

[Learn More](#)

Parameter	Value
Product Name	CAMK1G Rabbit pAb
Remarks/Alias	VWS1; CLICK3; CLICKIII; dj272L16.1; CAMK1G
Species	Human
Gene ID (Human)	57172
Gene ID	57172
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 217-476 of human CAMK1G (NP_065172.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96NX5
Protein Weight	53kDa
Shipping	Ice bag