


Anti-Camk2 Delta/Camk2 Gamma Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q13557 / Q13555

Item Number: CM0003270



Introduction

CamelBio offers an anti-CaMK2 delta/CAMK2 gamma rabbit monoclonal antibody, validated for WB, IP, and ELISA. It cross-reacts with human, mouse, and rat. The antibody targets a kinase crucial for cardiac Ca^{2+} homeostasis and excitation-contraction coupling.

[Learn More](#)

Parameter	Value
Product Name	CaMK2 delta/CAMK2 gamma Rabbit mAb
Remarks/Alias	CAM2; CAMK2; CAMKB; MRD54; CaMKIIbeta; CaMK2 delta/CAMK2 gamma
Species	Human
GeneID (Human)	817/818
GeneID	817/818
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human CaMK2 delta/CAMK2 gamma (Q13554)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13557/Q13555
Protein Weight	73 kDa
Shipping	Ice bag