


Anti-Camk2A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uqm7

Item Number: CM0020280



Introduction

High-quality rabbit polyclonal antibody targeting human CAMK2A (CaMKII α) for Western blot, IF/ICC, and ELISA. Validated in human, mouse, and rat. Essential for neuroscience research on synaptic plasticity and long-term potentiation.

[Learn More](#)

Parameter	Value
Product Name	CAMK2A Rabbit pAb
Remarks/Alias	CAMKA; MRD53; MRT63; CaMKIIalpha; CaMKIINalpha; CAMK2A
Species	Human
Gene ID (Human)	815
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human CAMK2A (NP_741960.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQM7
Protein Weight	54kDa
Shipping	Ice bag