

Anti-Camk4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16566

Item Number: CM0029331



Introduction

Rabbit polyclonal antibody targeting human CAMK4 (CaMKIV), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for neuroscience, immunology, and calcium signaling research.

[Learn More](#)

Parameter	Description
Product Name	CAMK4 Rabbit pAb
Remarks/Alias	caMK; CaMKIV; CaMK IV; CaMK-GR; CAMK4
Species	Human
Gene ID (Human)	814
Gene ID	814
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 304-473 of human CAMK4 (NP_001735.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16566
Protein Weight	52kDa
Shipping	Ice bag