

Anti-Camki Polyclonal Antibody For Wb, If/Icc, Elisa - Q14012

Item Number: CM0026803



Introduction

Rabbit polyclonal antibody targeting human CaMKI, suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse. Recognizes a 41kDa protein. Ideal for studying calcium signaling, neuronal differentiation, and cell cycle regulation.

[Learn More](#)

Parameter	Value
Product Name	CaMKI Rabbit pAb
Remarks/Alias	CAMKI; CaMKI
Species	Human
GeneID (Human)	8536
GeneID	8536
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 271-370 of human CaMKI (NP_003647.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14012
Protein Weight	41kDa
Shipping	Ice bag