


Anti-Camkii Delta (Camkd) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q13557

Item Number: CM0004372



Introduction

Rabbit monoclonal antibody targeting CaMKII delta (CAMKD) for Western blot, immunofluorescence, and ELISA. Cross-reacts with human and mouse. Suitable for studying cardiac calcium signaling and heart failure.

[Learn More](#)

Parameter	Value
Product Name	CaMKII delta Rabbit mAb
Remarks/Alias	CAMKD; CaMKII delta
Species	Human
GeneID (Human)	817
GeneID	817
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CaMKII delta (Q13557)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13557
Protein Weight	56 kDa
Shipping	Ice bag