


Anti-Camkii Delta Polyclonal Antibody For Wb, Elisa - Q13557

Item Number: CM0015392



Introduction

Rabbit polyclonal antibody against human CaMKII delta (CAMKD) for WB and ELISA applications. Cross-reacts with mouse. Useful for studying calcium signaling in cardiac and skeletal muscle.

[Learn More](#)

Parameter	Value
Product Name	CaMKII delta Rabbit pAb
Remarks/Alias	CAMKD; CaMKII delta
Species	Human
GeneID (Human)	817
GeneID	817
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human CaMKII delta (NP_742125.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q13557
Protein Weight	56kDa
Shipping	Ice bag