


Anti-Camkii Rabbit Monoclonal Antibody For Wb And Elisa - Q9Uqm7/Q13554/Q13557/Q13555

Item Number: CM0006696



Introduction

Rabbit monoclonal antibody targeting CaMKII (CAMK2A/B/D/G), validated for WB and ELISA. Recognizes human, mouse, and rat isoforms. Ideal for studying synaptic plasticity, calcium signaling, and kinase activity.

[Learn More](#)

Parameter	Value
Product Name	CaMKII Rabbit mAb
Remarks/Alias	CAMKA; MRD53; MRT63; CaMKIIalpha; CaMKIINalpha; CaMKII
Species	Human
Gene ID (Human)	815/816/817/818
Immunogen	A synthetic peptide corresponding to amino acids 200-300 of human CaMKII (NP_741960.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQM7/Q13554/Q13557/Q13555
Protein Weight	54kDa/72kDa/56kDa/62kDa
Shipping	Ice bag