


Anti-Camkii Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13554

Item Number: CM0003055



Introduction

Rabbit monoclonal antibody targeting CaMKII beta (CaMKII, CAMK2, CAMKB) (Q13554) for Western blot, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for neuroscience and muscle research.

[Learn More](#)

Parameter	Value
Product Name	CaMKII Rabbit mAb
Remarks/Alias	CAM2; CAMK2; CAMKB; MRD54; CaMKIibeta; CaMKII
Species	Human
GeneID (Human)	816
GeneID	816
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 550-650 of human CaMKII (Q13554)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13554
Protein Weight	73kDa
Shipping	Ice bag