


Anti-Camkii Rabbit Polyclonal Antibody For Wb And Elisa - Q13554

Item Number: CM0020344



Introduction

CaMKII Rabbit pAb validated for WB and ELISA, reacts with Human, Mouse, Rat. Targets CaMKII beta subunit involved in synaptic plasticity and calcium signaling. Recombinant immunogen, 73kDa. Suitable for neuroscience and muscle research.

[Learn More](#)

Parameter	Value
Product Name	CaMKII Rabbit pAb
Remarks/Alias	CAM2; CAMK2; CAMKB; MRD54; CaMKIIbeta; CaMKII
Species	Human
Gene ID (Human)	816
Gene ID	816
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13554
Protein Weight	73kDa
Shipping	Ice bag