


Anti-Camkk1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N5S9

Item Number: CM0023705



Introduction

Anti-CAMKK1 rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA in human and mouse samples. Ideal for studying calcium/calmodulin-dependent kinase signaling and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	CAMKK1 Rabbit pAb
Remarks/Alias	CAMKKA; CAMKK1
Species	Human
GeneID (Human)	84254
GeneID	84254
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 399-520 of human CAMKK1 (NP_757344.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N5S9
Protein Weight	56kDa
Shipping	Ice bag