


Anti-Camkk2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Rr4

Item Number: CM0029280



Introduction

Rabbit polyclonal antibody against human CAMKK2, validated for WB, IF/ICC, and ELISA. Detects human and mouse CAMKK2, a calcium/calmodulin-dependent kinase involved in neuronal signaling and AMPK activation. Suitable for neuroscience and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	CAMKK2 Rabbit pAb
Remarks/Alias	CAMKK; CAMKKB; CAMKK2
Species	Human
Gene ID	10645
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human CAMKK2 (NP_006540.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96RR4
Protein Weight	65kDa
Shipping	Ice bag