

Anti-Camkv Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Ncb2

Item Number: CM0023022

Introduction

Rabbit polyclonal antibody targeting human CAMKV (CaM kinase-like vesicle-associated protein) for Western blot (WB), immunofluorescence (IF/ICC), and ELISA. Detects human, mouse, and rat CAMKV.

[Learn More](#)

Parameter	Value
Product Name	CAMKV Rabbit pAb
Remarks/Alias	1G5; VACAMKL; CAMKV
Species	Human
Gene ID (Human)	79012
Gene ID	79012
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-285 of human CAMKV (NP_076951.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NCB2
Protein Weight	54kDa
Shipping	Ice bag