

Anti-Camp Polyclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P49913

Item Number: CM0023295



Introduction

Rabbit polyclonal antibody against human CAMP (cathelicidin antimicrobial peptide, LL37), validated for WB, IF/ICC, IF-P, IHC-P, ELISA. Reacts with human and mouse. For innate immunity and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	CAMP Rabbit pAb
Remarks/Alias	LL37; CAP18; CRAMP; HSD26; CAP-18; FALL39; FALL-39; CAMP
Species	Human
GeneID (Human)	820
GeneID	820
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 31-170 of human CAMP (NP_004336.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49913
Protein Weight	19kDa
Shipping	Ice bag