

Anti-Camp Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P49913

Item Number: CM0006169



Introduction

High-quality rabbit anti-CAMP monoclonal antibody for detection of human cathelicidin antimicrobial peptide (LL-37/hCAP-18). Validated in WB, IHC-P, and ELISA. Predicted MW 19kDa. Ideal for innate immunity and antimicrobial peptide research.

[Learn More](#)

Parameter	Value
Product Name	CAMP Rabbit mAb
Remarks/Alias	LL37; CAP18; CRAMP; HSD26; CAP-18; FALL39; FALL-39
Species	Human
GeneID (Human)	820
GeneID	820
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 71-170 of human CAMP (NP_004336.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P49913
Protein Weight	19kDa
Shipping	Ice bag