

Anti-Camp Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P51437

Item Number: CM0027493

Introduction

Rabbit polyclonal antibody against mouse CAMP (cathelicidin antimicrobial peptide, LL37). Validated for WB, IHC-P, and ELISA. Ideal for innate immunity and host defense research. UniProt P51437.

[Learn More](#)

Parameter	Value
Product Name	CAMP Rabbit pAb
Remarks/Alias	LL37; CAP18; CRAMP; HSD26; CAP-18; FALL39; FALL-39
Species	Mouse
Gene ID (Human)	820
Gene ID	12796
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 73-172 of mouse CAMP (NP_034051.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	P51437
Protein Weight	19 kDa
Shipping	Ice bag