


Anti-Cand1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q86Vp6

Item Number: CM0016049



Introduction

High-quality rabbit polyclonal antibody against human CAND1 (TIP120/TIP120A), validated for WB, IHC-P, and ELISA. Detects CAND1 in human, mouse, and rat samples. Suitable for ubiquitin-proteasome research.

[Learn More](#)

Parameter	Value
Product Name	CAND1 Rabbit pAb
Remarks/Alias	TIP120; TIP120A; CAND1
Species	Human
GeneID (Human)	55832
GeneID	55832
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 750-980 of human CAND1 (NP_060918.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86VP6
Protein Weight	136kDa
Shipping	Ice bag