

Anti-Cand2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O75155

Item Number: CM0019967



Introduction

High-quality rabbit polyclonal antibody against human CAND2, suitable for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying ubiquitin ligase complexes and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	CAND2 Rabbit pAb
Remarks/Alias	Tp120b; TIP120B; CAND2
Species	Human
Gene ID (Human)	23066
Gene ID	23066
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 850-1000 of human CAND2 (NP_001155971.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75155
Protein Weight	135kDa
Shipping	Ice bag