

Anti-Cant1 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Wvq1

Item Number: CM0025287



Introduction

High-quality rabbit polyclonal antibody against human CANT1 (soluble calcium-activated nucleotidase 1), validated for Western blotting (WB) and ELISA. Ideal for studying calcium-dependent nucleotide metabolism and proteoglycan synthesis in human samples.

[Learn More](#)

Parameter	Value
Product Name	CANT1 Rabbit pAb
Remarks/Alias	DBQD; EDM7; DBQD1; SCAN1; SHAPY; SCAN-1; CANT1
Species	Human
GeneID (Human)	124583
GeneID	124583
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 62-231 of human CANT1 (NP_620148.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8WVQ1
Protein Weight	45kDa
Shipping	Ice bag