


Anti-Cap1 Polyclonal Antibody For Wb, Elisa - Q01518

Item Number: CM0018125



Introduction

CAP1 Rabbit Polyclonal Antibody (Swiss-Prot Q01518) validated for Western blot (WB) and ELISA in human samples. Targets adenyl cyclase-associated protein 1, a key regulator of actin dynamics and cell polarity.

[Learn More](#)

Parameter	Value
Product Name	CAP1 Rabbit pAb
Remarks/Alias	CAP; CAP1-PEN; CAP1
Species	Human
Gene ID (Human)	10487
Gene ID	10487
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CAP1 (NP_006358.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q01518
Protein Weight	52kDa
Shipping	Ice bag