


Anti-Cap1 Rabbit Monoclonal Antibody For Wb, Elisa - Q01518

Item Number: CM0011164



Introduction

Anti-CAP1 Rabbit Monoclonal Antibody validated for WB and ELISA. Detects human CAP1 protein (~52 kDa). Supports research in actin dynamics, cell polarity, and mRNA localization.

[Learn More](#)

Parameter	Value
Product Name	CAP1 Rabbit mAb
Remarks/Alias	CAP; CAP1-PEN; CAP1
Species	Human
Gene ID (Human)	10487
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q01518
Protein Weight	52kDa
Shipping	Ice bag