


Anti-Capg Rabbit Polyclonal Antibody (Ko Validated) For Wb, Ihc-P, Elisa - P40121

Item Number: CM0030303



Introduction

KO validated rabbit polyclonal antibody against human macrophage-capping protein (CAPG). Suitable for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for actin dynamics and macrophage function studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CAPG Rabbit pAb
Remarks/Alias	MCP; AFCP; HEL-S-66; PG
Species	Human
Gene ID (Human)	822
Gene ID	822
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-348 of human CAPG (NP_001738.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40121
Protein Weight	38kDa
Shipping	Ice bag