

Anti-Capn10 Polyclonal Antibody For Wb, Elisa - Q9Hc96

Item Number: CM0028761

Introduction

CAPN10 Rabbit pAb is a polyclonal antibody raised against recombinant human CAPN10 protein. Suitable for WB and ELISA, it detects human calpain-10, a calcium-regulated protease implicated in glucose uptake and type 2 diabetes.

[Learn More](#)

Parameter	Value
Product Name	CAPN10 Rabbit pAb
Remarks/Alias	CANP10; NIDDM1; CAPN10
Species	Human
Gene ID (Human)	11132
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human CAPN10 (NP_075571.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HC96
Protein Weight	75kDa
Shipping	Ice bag