

Anti-Capn3 Polyclonal Antibody For Wb, If/Icc, Elisa - P20807

Item Number: CM0017813



Introduction

Rabbit polyclonal antibody targeting human CAPN3 (calpain-3/p94) for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for muscle proteolysis and limb-girdle muscular dystrophy research.

[Learn More](#)

Parameter	Value
Product Name	CAPN3 Rabbit pAb
Remarks/Alias	p94; CANP3; LGMD2; nCL-1; CANPL3; LGMD2A; LGMDD4; LGMDR1; CAPN3
Species	Human
GeneID (Human)	825
GeneID	825
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human CAPN3 (NP_775111.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P20807
Protein Weight	94kDa
Shipping	Ice bag