

Anti-Capn5 Polyclonal Antibody For Wb, If/Icc, Elisa - O15484

Item Number: CM0030584



Introduction

Rabbit polyclonal antibody against human CAPN5, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for studying calpain function and retinal diseases.

[Learn More](#)

Parameter	Value
Product Name	CAPN5 Rabbit pAb
Remarks/Alias	VRNI; ADNIV; HTRA3; nCL-3; CAPN5
Species	Human
GeneID (Human)	726
GeneID	726
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 361-640 of human CAPN5 (NP_004046.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15484
Protein Weight	73kDa
Shipping	Ice bag