


Anti-Capn7 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y6W3

Item Number: CM0019784



Introduction

High-quality rabbit polyclonal antibody against human CAPN7, suitable for WB, IF/ICC, and ELISA. Validated in human, mouse, and rat. Ideal for studying calcium-regulated thiol-protease in the nucleus.

[Learn More](#)

Parameter	Value
Product Name	CAPN7 Rabbit pAb
Remarks/Alias	PALBH; CALPAIN7; CAPN7
Species	Human
GeneID (Human)	23473
GeneID	23473
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-220 of human CAPN7 (NP_055111.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6W3
Protein Weight	93kDa
Shipping	Ice bag