

Anti-Capn9 Polyclonal Antibody For Ihc-P, Elisa - O14815

Item Number: CM0021883



Introduction

Rabbit polyclonal antibody against human CAPN9 (calpain-9, nCL-4), validated for immunohistochemistry (paraffin) and ELISA. Ideal for stomach-specific protease research.

[Learn More](#)

Parameter	Value
Product Name	CAPN9 Rabbit pAb
Remarks/Alias	GC36; nCL-4; CAPN9
Species	Human
Gene ID (Human)	10753
Gene ID	10753
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 391-690 of human CAPN9 (NP_006606.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	O14815
Protein Weight	79kDa
Shipping	Ice bag