

Anti-Capn9 Rabbit Polyclonal Antibody For Wb, Elisa - O14815

Item Number: CM0014714



Introduction

Anti-CAPN9 rabbit polyclonal antibody for WB, ELISA. Recognizes human, mouse, rat CAPN9, a digestive tract-specific calpain (nCL-4) involved in calcium-regulated proteolysis.

[Learn More](#)

Parameter	Value
Product Name	CAPN9 Rabbit pAb
Remarks/Alias	GC36; nCL-4; CAPN9
Species	Human
Gene ID (Human)	10753
Gene ID	10753
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 391-690 of human CAPN9 (NP_006606.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14815
Protein Weight	79kDa
Shipping	Ice bag