


Anti-Caprin2 Rabbit Polyclonal Antibody For Wb, Elisa - Q6Imn6

Item Number: CM0022291



Introduction

High-quality anti-CAPRIN2 rabbit polyclonal antibody validated for WB and ELISA. Detects human and rat CAPRIN2 (EEG1, C1QDC1), a key regulator of Wnt signaling and RNA metabolism. Ideal for cancer, neurobiology, and cell signaling research.

[Learn More](#)

Parameter	Specification
Product Name	CAPRIN2 Rabbit pAb
Remarks/Alias	EEG1; EEG-1; C1QDC1; RNG140; CAPRIN2
Species	Human
GeneID (human)	65981
GeneID	65981
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q6IMN6
Protein Weight	126kDa
Shipping	Ice bag