

Anti-Capza1 Rabbit Polyclonal Antibody For Wb, Elisa - P52907

Item Number: CM0026630



Introduction

CAPZA1 Rabbit polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, rat CAPZA1, the F-actin-capping protein subunit alpha-1, involved in actin filament dynamics and cell junctions.

[Learn More](#)

Parameter	Value
Product Name	CAPZA1 Rabbit pAb
Remarks/Alias	CAPZ; CAZ1; CAPP1; CAPZA1
Species	Human
Gene ID (human)	829
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-286 of human CAPZA1 (NP_006126.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52907
Protein Weight	33kDa
Shipping	Ice bag