

Anti-Capza2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P47755

Item Number: CM0021778

Introduction

CAPZA2 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Detects human and mouse CAPZA2, an F-actin-capping protein. Suitable for actin cytoskeleton and dynamics research.

[Learn More](#)

Parameter	Value
Product Name	CAPZA2 Rabbit pAb
Remarks/Alias	CAPZ; CAPPA2; CAPZA2
Species	Human
GeneID (Human)	830
GeneID	830
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-286 of human CAPZA2 (NP_006127.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P47755
Protein Weight	33kDa
Shipping	Ice bag