

Anti-Capzb Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P47756

Item Number: CM0027981



Introduction

Rabbit polyclonal antibody targeting CAPZB, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat F-actin-capping protein subunit beta (31 kDa). Ideal for cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	CAPZB Rabbit pAb
Remarks/Alias	CAPB; CAPZ; CAPPB; CAPZB
Species	Human
GeneID	832 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-272 of human CAPZB (NP_004921.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P47756
Protein Weight	31kDa
Shipping	Ice bag