


## Anti-Carbonic Anhydrase 1 (Ca1) Rabbit Polyclonal Antibody For Wb And Elisa - P00915

Item Number: CM0014881



### Introduction

Rabbit polyclonal antibody targeting human Carbonic Anhydrase 1 (CA1), validated for Western blotting and ELISA. Cross-reacts with human and mouse. Designed for research into carbon dioxide hydration, pH homeostasis, and related metabolic pathways.

[Learn More](#)

| Parameter        | Value                                                                                                                              |
|------------------|------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Carbonic Anhydrase 1 (CA1) Rabbit pAb                                                                                              |
| Remarks/Alias    | CAB; CA-I; Car1; HEL-S-11; Carbonic Anhydrase 1 (CA1)                                                                              |
| Species          | Human                                                                                                                              |
| GeneID           | 759                                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-261 of human Carbonic Anhydrase 1 (NP_001122303.1) |
| Source           | Rabbit                                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                                              |
| Application      | WB, ELISA                                                                                                                          |
| Cross Reactivity | Human, Mouse                                                                                                                       |
| SWISS            | P00915                                                                                                                             |
| Protein Weight   | 29kDa                                                                                                                              |
| Shipping         | Ice bag                                                                                                                            |