


Anti-Carbonic Anhydrase 9 (Ca9/G250) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q16790

Item Number: CM0004582



Introduction

Anti-Carbonic Anhydrase 9 (CA9/G250) Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Specific to human CA9, a tumor-associated antigen. Ideal for cancer research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Carbonic Anhydrase 9 (CA9/G250) Rabbit mAb
Remarks/Alias	MN; CAIX; Carbonic Anhydrase 9 (CA9/G250)
Species	Human
GeneID (Human)	768
GeneID	768
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Carbonic Anhydrase 9 (CA9/G250) (Q16790)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q16790
Protein Weight	50kDa
Shipping	Ice bag