


Anti-Carbonic Anhydrase 9 (CA9/G250) Rabbit Polyclonal Antibody For Wb - Q16790

Item Number: CM0013952



Introduction

Rabbit polyclonal antibody specific to human Carbonic Anhydrase 9 (CA9/G250, Q16790). Suitable for Western blot, IHC-P, IF/ICC, and ELISA. Recombinant protein immunogen. Useful for researching tumor hypoxia and pH regulation.

[Learn More](#)

Parameter	Value
Product Name	Carbonic Anhydrase 9 (CA9/G250) Rabbit pAb
Remarks/Alias	MN; CAIX; Carbonic Anhydrase 9 (CA9/G250)
Species	Human
GeneID	768 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-151 of human Carbonic Anhydrase 9 (CA9/G250) (NP_001207.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q16790
Protein Weight	50 kDa
Shipping	Ice bag