

Anti-Carboxypeptidase M (Cpm) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P14384

Item Number: CM0026996



Introduction

High-quality rabbit polyclonal antibody against human carboxypeptidase M (CPM). Validated for WB, IHC-P, and ELISA. Suitable for studying peptide hormone regulation and membrane-localized proteolysis. Cross-reacts with human and rat.

[Learn More](#)

Parameter	Value
Product Name	CPM Rabbit pAb
Remarks/Alias	CPM
Species	Human
GeneID (Human)	1368
GeneID	1368
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 204-443 of human CPM (NP_001005502.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P14384
Protein Weight	51kDa
Shipping	Ice bag