


Anti-Car/Cxadr Polyclonal Antibody For Wb, Elisa - P78310

Item Number: CM0019153



Introduction

Rabbit polyclonal antibody against human CAR/CXADR, validated for WB and ELISA; cross-reacts with mouse and rat. Ideal for studying tight junctions, viral receptors, and leukocyte migration.

[Learn More](#)

Parameter	Value
Product Name	CAR/CXADR Rabbit pAb
Remarks/Alias	CAR; HCAR; CAR4/6; CAR/CXADR
Species	Human
GeneID (Human)	1525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-237 of human CAR/CXADR (NP_001329.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78310
Protein Weight	40kDa
Shipping	Ice bag