


Anti-Car/Cxadr Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P78310

Item Number: CM0007253



Introduction

Recombinant rabbit monoclonal antibody against human CAR/CXADR, a tight junction and viral receptor. Validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studies in cell adhesion, epithelial biology, and host-virus interactions.

[Learn More](#)

Parameter	Value
Product Name	CAR/CXADR Rabbit mAb
Remarks/Alias	CAR; HCAR; CAR4/6; CAR/CXADR
Species	Human
Gene ID (Human)	1525
Gene ID	1525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-237 of human CAR/CXADR (NP_001329.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78310
Protein Weight	40kDa
Shipping	Ice bag