

Anti-Card10 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bwt7

Item Number: CM0030634

Introduction

Rabbit polyclonal antibody against CARD10 (CARMA3, BIMP1), validated for western blot, IHC-P, and ELISA. Detects human, mouse, and rat CARD10. Suitable for NF-kB pathway studies.

[Learn More](#)

Parameter	Value
Product Name	CARD10 Rabbit pAb
Remarks/Alias	BIMP1; IMD89; CARMA3; CARD10
Species	Human
GeneID (Human)	29775
GeneID	29775
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 863-1032 of human CARD10 (NP_055365.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BWT7
Protein Weight	116kDa
Shipping	Ice bag