


Anti-Card11 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bxl7

Item Number: CM0025565



Introduction

High-quality rabbit polyclonal antibody against CARD11 (CARMA1, Caspase recruitment domain-containing protein 11), suitable for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for adaptive immunity and NF-kB signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CARD11 Rabbit pAb
Remarks/Alias	PPBL; BENTA; BIMP3; IMD11; CARMA1; IMD11A; CARD11
Species	Human
GeneID	84433 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 985-1154 of human CARD11 (NP_115791.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXL7
Protein Weight	133kDa
Shipping	Ice bag